

Kinetic

 WATCH NOW

Scan me on your
phone camera
to watch our
assembly video.



USER MANUAL

PRODUCT CODE: WC02010019

INTRODUCING YOUR KINETIC COMPACT WHEELCHAIR

Foldable and effortlessly transportable, the Kinetic Compact Wheelchair is the ideal companion designed to provide you with maximum comfort whilst increasing independence in your daily life.

Due to CareCo's commitment to continuous product improvement, product information and designs are subject to change without notice. No claims can be made based on the data or illustrations in this user manual.



CareCo (UK) Limited, 1 Turing Court,
Great Notley, Braintree, CM77 7AT

Tel: 0333 015 5000
www.CareCo.co.uk

Thank you for choosing CareCo.



This manual can be viewed in PDF format at the following URL: **www.careco.co.uk**

BEFORE USE

Before using your wheelchair, please ensure you, and any persons that may assist you in its use, take the time to read, understand and follow the instructions given within this user manual.

We recommend that you keep this user manual in a safe place for future reference

If you have any questions about product use, maintenance, or safety, please contact CareCo.

KEY DEFINITIONS

This user manual contains warnings and suggestions marked with the following symbols:



WARNING! Highlights a potentially hazardous condition/situation for the user and/or product.



SUGGESTION General suggestions and advice for operating the product.

Need help with assembly?

We've put together a step-by-step assembly guide to help you get your product set up safely and correctly. Just scan the QR code with your phone camera to play the video.



CONTENTS

GETTING STARTED	6
CARTON CONTENTS	
PRODUCT COMPONENTS	
SETTING UP YOUR KINETIC COMPACT WHEELCHAIR	
HOW TO USE YOUR WHEELCHAIR	9
USER INSTRUCTIONS	
APPLYING THE BRAKES	
LEGRESTS	
FOOTPLATES	
SEAT & BACK COVERS	
SEAT BELT	
STORAGE SPACE	
ADDITIONAL STORAGE SPACE	
FOLDING, TRANSPORTING & STORING YOUR WHEELCHAIR	15
FOLDING INSTRUCTIONS	
TRANSPORTATION	
STORING THE WHEELCHAIR	
COMFORT ADJUSTMENTS	18
FOOTPLATE HEIGHT ADJUSTMENT	
LOOKING AFTER YOUR WHEELCHAIR	19
CLEANING	
MAINTENANCE	
MAINTENANCE SCHEDULE	
REPAIRS	

OTHER ADJUSTMENTS	20
BRAKE ADJUSTMENT	
WARNINGS & TECHNICAL INFORMATION	21
SPECIFICATIONS	
PRODUCT REISSUING	
RECYCLING & DISPOSAL	
MATERIALS	
INTENDED USE	
INTENDED USER	
INTENDED ENVIRONMENT	
PRODUCT IDENTIFICATION	
GENERAL WARNINGS	
EXPECTED LIFETIME	25
WARRANTY	26
NOTES	27

GETTING STARTED

ITEM NUMBER	DESCRIPTION	QUANTITY
1	Kinetic Compact Wheelchair	1
2	Legrest Set	1
3	Carry Case	1
4	User Manual	1
5	Hex Spanner	1

Contact CareCo immediately if any of the listed contents are missing.

We strongly recommend that you keep the box and internal packaging for the duration of your warranty period.

PRODUCT COMPONENTS



- 1. Push Handle
- 2. Half-Folding Back Mechanism
- 3. Attendant Brake
- 4. Armrest
- 5. Lap Belt
- 6. Leg Rest
- 7. Rear Wheel
- 8. Front Castor

SETTING UP YOUR KINETIC COMPACT WHEELCHAIR



1. Remove the carry case from the outer packaging.

2. Remove all wheelchair components from the carry case and place them onto a level surface.



3. Place the wheelchair upright as shown.



4. Push downwards on the corners of the wheelchair seat to open the wheelchair.



5. Lift the backrest handles upwards.

An audible “click” will confirm the handles are secure.



6. Lift the armrests, then locate the spring pins into the lower wheelchair frame.



7. Place the velcro seat and back covers onto the wheelchair.




8. Place the legrests onto the wheelchair frame at a 90-degree angle, then rotate inwards to secure.

HOW TO USE YOUR WHEELCHAIR

USER INSTRUCTIONS

TRANSFERRING TO & FROM

Always ensure the brakes are applied before getting in/out of the wheelchair.




TRANSFER TO

Position your body as close to the wheelchair as possible. Support your body using the armrests, before lowering into the seat.

TRANSFER FROM

Support your body using the armrests and place your feet on the ground. Push upwards to lift out of the seat.

GRADIENTS



ASCENDING

Tilt your upper body forward when travelling up gradients for increased stability.




DESCENDING

Tilt your upper body back when travelling down gradients for increased stability.


THRESHOLDS & FIXED OBSTACLES

Approach thresholds vertically to avoid wheelchair misbalance.

Incorrect Approach



Correct Approach



SUGGESTION Use caution when operating the wheelchair near kerbs, ledges, and doorway thresholds.

DOORS & GATES

Determine whether the door/gate opens towards or away from the user.

Reach for the handle, then carefully move forwards or backwards depending on the door/gate orientation.

STAIRS, ESCALATORS AND LIFTS

The wheelchair is not designed to travel up/down stairs and escalators. Use a lift whenever possible.

LIFTING & CARRYING

Lift/carry the wheelchair from areas of the frame that are clear from pinch points.

APPLYING THE BRAKES

ATTENDANT BRAKE

STOPPING

Gradually squeeze both brake levers to stop the wheelchair. Release the brake levers to start moving again.



TURNING

Squeeze the left and right brake levers individually to rotate the wheelchair.

Release the brake lever to stop rotating.


LOCKING


Locate the locking lever in the handle catch to lock the attendant brake.

UNLOCKING

Lift the locking lever to unlock the attendant brake.



 **WARNING!** Ensure the wheelchair is stationary when locking the brakes.

 **SUGGESTION** Regularly test the brakes to ensure they function correctly.

LEGRESTS



ATTACHING

Place the legrest onto the frame at a 90-degree angle, then rotate inwards.

An audible “click” will confirm the legrest is secure.

Repeat on the other side.



REMOVING

Lift the legrest lever to release the leg rest from the frame, then rotate outwards.

Lift the legrest away from the frame.

Repeat on the other side.

FOOTPLATES

LIFTING

Fold the footplates upwards to assist the user when leaving the wheelchair and when folding the wheelchair.

LOWERING

Fold the footplates downwards to aid the user when travelling in the wheelchair.



SEAT & BACK COVERS

Your Kinetic Compact Wheelchair seat and backrest covers are removeable and can be washed if required.

ATTACHING

Place the covers onto the wheelchair back & seat canvas.



REMOVING

Pull the covers away from the wheelchair back & seat canvas.


SEAT BELT

FASTENING

Insert the male buckle into the female buckle.

RELEASING

Press the exposed sides of the male buckle while pulling away from the centre



SUGGESTION Regularly test and adjust the seat belt to suit the user.

STORAGE SPACE

Your Kinetic Compact Wheelchair features a rear storage pocket, ideal for holding personal belongings and small accessories.

MAX CARRY WEIGHT: 5KG



ADDITIONAL STORAGE SPACE



CARECO WHEELCHAIR BAG

Attach additional storage space to the back of your wheelchair by looping the bag straps over the wheelchair's push handles.



POPEM PANNIER BAG

Attach additional storage space to the side of your wheelchair by looping the bag straps over the wheelchair's armrest.

FOLDING, TRANSPORTING & STORING YOUR WHEELCHAIR

FOLDING INSTRUCTIONS



1. Remove the velcro seat & backrest covers from the wheelchair.



2. Pull the legrest lever(s) and rotate the legrest(s) outwards.

Lift the legrest(s) upwards to remove from the frame.



3. Pull the armrest spring pins, then lower the armrests.



4. Press the backrest levers, then fold down the backrest handles.



5. Bring the corners of the wheelchair seat and backrest inwards to fold the wheelchair.



6. Place the wheelchair components into the carry case.



7. Zip the carry case, ready for transport or storage.



WARNING! Keep hands and body from the cross-frame pivot to avoid possible pinch points.

TRANSPORTATION

STOWING THE WHEELCHAIR

Ensure the brakes have been applied and any accessories have been secured before transporting.

Fold the wheelchair, then insert it into the carry bag.

Position the bag appropriately within the vehicle to prevent excess movement during transit.



WARNING! This product is not suitable as a seat in a motor vehicle.



TRANSPORTING BY AIR & RAIL

Prior to your journey, contact the service provider to confirm their transport-specific requirements.

STORING THE WHEELCHAIR

SHORT TERM STORAGE

Store the product in a cool and dry environment.

LONG TERM STORAGE

Use the maintenance schedule found within this user manual as a guide before and after long term storage.

Store the product in a cool and dry environment.



WARNING! Storing the product in extreme temperatures can cause premature rusting & damage.

COMFORT ADJUSTMENTS

FOOTPLATE HEIGHT ADJUSTMENT

Difficulty Level

LOW

Moderate

High



Tools Required

Hex Spanner (included)

1. Using a spanner, loosen the hex head at the base of the footrest.
2. Lift or lower the footplate until it is at the desired height.
3. Tighten the hex head to secure the position.
Repeat on the other side.



LOOKING AFTER YOUR WHEELCHAIR

CLEANING

Regular use of the wheelchair will accumulate dirt and debris over time.

Clean the wheelchair regularly using mild soap and warm water, then dry using a cloth.

 **WARNING!** Do not use bleach or strong chemicals to clean the wheelchair.

 **SUGGESTION** Regular cleaning & maintenance can help prolong the wheelchairs lifespan.

MAINTENANCE

The wheelchair is made up of components which over time will become susceptible to wear and tear. Knowing how to properly care and maintain your product will not only ensure peak performance and safe operation but will also help maintain operation for years to come.

MAINTENANCE SCHEDULE

Carry out regular inspections to keep the wheelchair in optimum condition.

COMPONENT	INSPECTION DESCRIPTION	AT ANY TIME	WEEKLY	MONTHLY	QUARTERLY
FASTENERS	Ensure the fasteners are tight and components are not loose.				○
WHEELS & TYRES	Ensure tyres have sufficient tread and the wheels are in a good condition.				○
BRAKES	Ensure the brakes function as intended.	○			
UPHOLSTERY	Check for wear and tear.		○		
OVERALL	Overall check for wear and residue build-up. Cleaning if required.			○	

REPAIRS

Certain wheelchair components can be repaired by the user with guidance provided by CareCo. Other components can only be repaired/replaced by a service technician.

If for any reason the wheelchair requires repair due to a defect, contact CareCo immediately.

OTHER ADJUSTMENTS

BRAKE ADJUSTMENT

Difficulty Level

LOW

Moderate

High



Tools Required

8mm Spanner x1

10mm Spanner x1

The user can adjust the brake using simple household tools.

Loosen the lock nut (B) by rotating it anticlockwise.

Adjust the brake cable tension by rotating the brake adjuster (A).

Tighten the lock nut (B) by rotating it clockwise.



WARNINGS & TECHNICAL INFORMATION

SPECIFICATIONS

OVERALL LENGTH (MM)	900
OVERALL WIDTH (MM)	600
OVERALL HEIGHT (MM)	930
FOLDED LENGTH (MM)	330
FOLDED WIDTH (MM)	310
FOLDED HEIGHT (MM)	740
FRONT WHEEL DIAMETER (IN)	8
REAR WHEEL DIAMETER (IN)	12
CARRY WEIGHT (KG)	9.5
WEIGHT OF HEAVIEST COMPONENT (KG)	4.5
WEIGHT CAPACITY (KG)	100
SEAT WIDTH (MM)	440
SEAT DEPTH (MM)	430
SEAT HEIGHT (MM)	480
MAX/MIN OPERATING TEMPERATURE (°C)	40°C / -5°C
MAX/MIN STORAGE TEMPERATURE (°C)	40°C / -5°C
TYRE TYPE	FRONT: PVC – REAR: PU
SUITABLE AS A SEAT IN A MOTOR VEHICLE (YES/NO)	NO

PRODUCT REISSUING

Ensure the product is safe and functions as intended before re-issuing.

RECYCLING & DISPOSAL

Recycling or disposal of the wheelchair should be completed through an approved distributor or disposal facility. (Local disposal & recycling regulations may apply)

Alternatively, contact CareCo for advice on disposal.

MATERIALS:

WHEELCHAIR PACKAGING:

Cardboard Box, PE high density polyethylene bag.

WHEELCHAIR UPHOLSTERY:

Nylon Fabric, Foam.

PLASTIC:

Hand grip, attendant brake handle, front castor, rear wheel, footplate, armrest pad.

ALUMINIUM:

Frame, handle tube, footrest tube, cross brace, front castor fork.

STEEL:

Bearing, fasteners, rear wheel brake.

INTENDED USE

This product is intended for users who experience difficulty in walking and/or have reduced mobility.

INTENDED USER


To aid a user with mobility limitations to move around safely.

INTENDED ENVIRONMENT

This product is designed to be used both indoors and outdoors on firm, level ground. Take care when using the product in inclement weather conditions.

PRODUCT IDENTIFICATION

The product identification label is located on the cross frame.




DEFINITIONS

REF	Product code		Electronic manual
MD	Medical device		Operator's manual
SN	Serial number		Manufacturer's address
UDI	Unique device identification		UKCA mark
	Date of manufacturer		

GENERAL WARNINGS

Read the following carefully to avoid damage and/or injury caused by improper operation.

- Do not exceed the maximum user weight indicated on the frame label.
- Do not stand on the product or the footrests at any time.
- Place the product on a stable, level surface before getting in and out.
- Ensure the products castors face forward before getting in and out.
- Apply the brakes when using the product in an elevator, wheelchair lift and when getting in and out.
- When ascending/descending stairs, ensure the attendant holds on to secured parts only.
- Potholes and unevenness in the ground can cause the product to tip over.
- Ensure tyres have sufficient tread.
- If the product is exposed to direct sunlight/ cold temperature for a long time, the components can become hot or cold.
- Avoid travelling on roads and use pavements whenever possible.
- Do not use the product if any of the parts are not functioning correctly or are damaged.
- A sudden change of direction may cause the product to tip over
- Do not attach anything to the wheels.
- Periodically test the brakes to ensure safe use.
- Beware of pinch points when using the wheelchair.
- This product cannot be used as a seat in a motor vehicle. 

EXPECTED LIFETIME

The Expected Lifetime is the estimated period of time for which the product is expected to function according to its intended use.

Various items such as weight capacity, user behaviour, external conditions and intensity of use can affect the expected lifetime.

WEIGHT CAPACITY:

If the user exceeds the weight limit specified, it may affect the product's expected lifetime.

USER BEHAVIOUR:

Using excessive force, rough handling or dropping components can significantly shorten the products's expected lifetime.

EXTERNAL CONDITIONS:

Extreme temperatures, humidity and exposure to harsh weather conditions can accelerate the deterioration of the product.

USAGE:

A product that is used daily may experience more wear and tear compared to a product that is used occasionally.

The expected lifetime for the wheelchair is 5 years providing the instructions and information within this user manual are followed.

Regular cleaning, lubrication of moving parts, and tightening of loose fixings will prevent premature wear and damage and extend the product's lifetime.

The expected lifetime is not related to the product warranty.

WARRANTY

If a component requires repair or replacement due to a specific manufacturing or material defect within its warranty period, it will be repaired or replaced free of charge by CareCo.

- **FRAME: 1-YEAR LIMITED WARRANTY**
- **DURABLE PARTS: 1-YEAR LIMITED WARRANTY**

Consumables are not covered under normal warranty circumstances unless they have clearly suffered excessive wear or damage as a direct result of a manufacturing defect. These items include grips, tyres, and rubber accessories.

The guarantee for the wheelchair is not transferable.

Parts that have been previously replaced that require subsequent repair or replacement due to manufacturing defects will be honoured for the remainder of the wheelchair's warranty period.

Parts that have been purchased by the end user after the original warranty period has expired have a warranty period of 12 months.

CareCo is not responsible if the wheelchair requires repair or replacement as a direct result of:

- | | |
|---|---|
| • Changes or modifications made to the product that fall outside of the manufacturer's specifications. | • Negligence, Accidental Damage, and Improper use. |
| • The product not being serviced or maintained according to the manufacturer's recommendations contained in this user manual. | • Repairs attempted before CareCo have been notified. |
| | • The product being purposely damaged or mistreated. |

In the unlikely event that this product is recalled, CareCo will notify the user.

Please ensure that your contact information is up to date.

Alternatively, please visit our website for any product alerts: www.CareCo.co.uk



CareCo (UK) Limited, 1 Turing Court,
Great Notley, Braintree, CM77 7AT, England



SELF PROPELLED
WHEELCHAIR



24 05 02