

Directions for use

Step 1

Align the middle of your foot with the back legs of the walker.

Step 2

Lift the walker up and place it comfortably in front of you, making sure all four legs end up on the floor.

Step 3

Carefully, walk forward one step at a time.

Folding frame



Wheelchair to Walker

Making certain the wheels on the wheelchair are locked lift your body by pushing down on the armrest and carefully transfer one hand at a time to the walker.

Custom Adjustments

With your arm extended down, the handle should be at wrist height. With your hand placed on the hand-grip, your elbow should be flexed 20° - 30°.

Caution:

Be aware of the following:

- Electrical cords in the way
- Liquids on the floor
- Throw rugs
- Toys or debris
- Furniture
- Doors and door steps
- Gaps in floor

Warning:

1. Use wheels on front legs only or on all four legs with appropriate rear wheel brakes.
2. Do not push one side of the walker ahead and then other side without lifting the walker first.
3. Use Drive accessories only. Use of other manufacturers' accessories can result in improper fit and risk fall and/or injury.
4. Do not use on stairs or escalators.