

CONNECT

ZORA

*The world's
lightest
powerchair*



Li
LITHIUM
BATTERY

Crafted from carbon fibre for optimum strength and durability.

- Weighs only 13kg
- Airline-friendly batteries
- Full carbon fibre frame
- Innovative control panel
- Quick and easy folding
- Magnesium alloy wheels



YEAR FRAME WARRANTY



Expert engineering and the latest technology means this powerchair weighs less than any other you'll find on the market.

Five-times stronger than steel and twice as stiff, it also resists corrosion and it doesn't fatigue over time. The higher compression strength of carbon fibre coupled with the reinforced seat panel means that despite its lightweight build, this powerchair still supports up to 120 kg in weight. In short? The Zora is the way to go if you want a powerchair that will last.

Dimensions Height x Length x Width	86 x 90 x 58cm
Folded Dimensions Height x Length x Width	68.5 x 72 x 30cm
Range 1 battery	Up to 15 km
Range 2 batteries	Up to 30 km
Product Weight	13 kg
Max User Weight	120 kg