

ErgoPro

Truly a game changer for ergonomic wheelchairs.



**Exceptional
build quality**



Produced to the highest standards, with a 5-year warranty for peace of mind



**A smooth
mover**



Solid tyres reduce rolling resistance while the anti-tip wheels give you confidence



**Customised for
your comfort**



Tension-adjustable backrest for optimum support and pressure relief

ErgoPro

The CareCo Ergo Pro Wheelchairs are more than just your everyday mobility aid. They're modern, sleek and slick, with quick-release wheels, a tension adjustable backrest and a two-tone metallic finish, as well as a 5-year warranty on the frame for total peace of mind while you're on the move. Puncture-proof tyres offer minimal rolling resistance for a smooth ride, anti-tip wheels boost its safety credentials on kerbs and hills.




- Two-tone metallic finish
- Puncture-proof tyres
- Easily folds
- Removable footrests
- Ergonomic design




A frame to suit your frame

Transit



Seat Width	16"	18"	20"	22"
Height (mm)	935	935	935	935
Length (mm)	1090	1090	1090	1090
Width (mm)	575	625	670	740
Folded Height (mm)	745	765	765	850
Folded Length (mm)	1090	1090	1090	1090
Folded Width (mm)	315	315	315	350
Carry Weight (kg)	11.4	12.1	12.8	13.5
Maximum User Weight	115Kg	115Kg	135Kg	160kg

Self Propelled



Seat Width	16"	18"	20"	22"
Height (mm)	935	935	935	935
Length (mm)	1110	1110	1110	1110
Width (mm)	563	625	687	749
Folded Height (mm)	745	765	765	850
Folded Length (mm)	1110	1110	1110	1110
Folded Width (mm)	325	295	315	315
Carry Weight (kg)	10	11.14	11.85	12.65
Maximum User Weight	115Kg	115Kg	135Kg	160kg