



NEW



MINIMUS

MICROMUS

MAXIMUS

RANGE
8 miles

TOTAL WEIGHT
19.4 kg

CARRY WEIGHT
17.8 kg

TILLER BAR
Straight

TILLER MATERIAL
Aluminium

WHEEL SIZE
6.5"

DIMENSIONS
height **71-91.0 cm**
length **90.5 cm**
width **42.0 cm**

SEAT WIDTH
13"

OFF BOARD CHARGING
Yes

MAX USER WEIGHT
115 kg

SWIVEL SEAT
No

FRONT & REAR LIGHTS
No

USB CHARGING PORT
No

IDEAL TERRAIN
Pavement

COLOURS

Lime

Azure

CARRY CASE
Optional Extra

ARMREST
Optional Extra

BASKET
Optional Extra

AIRLINE FRIENDLY
Yes

RANGE
8 miles

TOTAL WEIGHT
17.0 kg

CARRY WEIGHT
15.4 kg

TILLER BAR
Straight

TILLER MATERIAL
Carbon Fibre

WHEEL SIZE
5.8"

DIMENSIONS
height **86.0 cm**
length **90.0 cm**
width **43.0 cm**

SEAT WIDTH
16"

OFF BOARD CHARGING
Yes

MAX USER WEIGHT
100 kg

SWIVEL SEAT
No

FRONT & REAR LIGHTS
No

USB CHARGING PORT
No

IDEAL TERRAIN
Pavement

COLOURS

Indigo

Dune

CARRY CASE
No

ARMREST
No

BASKET
Optional Extra

AIRLINE FRIENDLY
Yes

RANGE
24 miles

TOTAL WEIGHT
29.9 kg

CARRY WEIGHT
21.5 kg

TILLER BAR
Delta

TILLER MATERIAL
Aluminium

WHEEL SIZE
7.5"

DIMENSIONS
height **84.0 cm**
length **95.6 cm**
width **49.0 cm**

SEAT WIDTH
16"

OFF BOARD CHARGING
Yes

MAX USER WEIGHT
115 kg

SWIVEL SEAT
Yes

FRONT & REAR LIGHTS
Yes

USB CHARGING PORT
Yes

IDEAL TERRAIN
Pavement & Bumpier ground

COLOURS

Crimson

Onyx

CARRY CASE
No

ARMREST
Included

BASKET
Included

AIRLINE FRIENDLY
Optional Extra

VS

VS