

Finestra



kondor



Finestra VS kondor

COLOURWAY	Bronze	Silver/Blue
MAX SPEED	4MPH/8MPH SWITCHABLE	4MPH/6MPH SWITCHABLE
BATTERY TYPE	LITHIUM-IRON PHOSPHATE	LEAD-ACID
CHARGING TIME	6-8 HOURS	12-16 HOURS
MAXIMUM POTENTIAL RANGE	75KM	44KM
TYRE TYPE	PNEUMATIC	PNEUMATIC
WHEEL DIAMETER	15IN	15IN
TOTAL WEIGHT	181KG	229KG
WEIGHT CAPACITY	200KG	200KG
EXTERIOR CHARGING PORT	✓	
USB A & C CHARGING PORT	✓	
INTERIOR CHARGING PORT		✓
EXTERIOR STORAGE COMPARTMENT	✓	
INTERIOR STORAGE COMPARTMENT	✓	✓
WINDSCREEN WIPER	✓	✓
CABIN HEATER	✓	✓
CABIN LIGHT	✓	✓
REVERSING CAMERA	✓	✓
LOCKABLE CABIN DOORS	✓	✓
FRONT & REAR LIGHTING	✓	✓
WING MIRRORS	✓	✓
ADJUSTABLE CAPTAIN SEAT	✓	✓
FRONT & REAR SUSPENSION	✓	✓