

ZORA



ZR1



ZORA EXPLORER



ZORA

ZR1

ZORA

EXPLORER

Lithium-Ion  
(Includes 1x Additional Battery)

BATTERY TYPE	Lithium-Ion	Lithium-Ion	Lithium-Ion (Includes 1x Additional Battery)
JOYSTICK TYPE	Analog	Digital LCD	Analog
CHARGING TIME	6-8 Hours	6-8 Hours	6-8 Hours
MAXIMUM POTENTIAL RANGE (KM)	15	16.4	30
TYRE TYPE	Solid	Solid	Solid
TOTAL WEIGHT	13	14	18
WEIGHT CAPACITY	120	150	136
OFF-BOARD CHARGING	✓	✓	✓
SIDE STORAGE BAG	✓	✓	✓
FOLDABLE ARMRESTS	✓	✓	✓
SAFETY BELT	✓	✓	✓
ANTI-TIP WHEELS	✓	✓	✓
SUSPENSION			✓
FRONT WHEEL DIAMETER (IN)	6.5	6.5	7
REAR WHEEL DIAMETER (IN)	10.6	11	12
LXWXH (MM)	900x580x860	990x600x980	900x603x875
FOLDED LXWXH (MM)	720x295x685	720x300x720	710x318x690
SEAT HEIGHT (MM)	460	515	480
SEAT WIDTH (MM)	420	450	420
SEAT DEPTH (MM)	400	395	400
MATERIAL	Carbon Fibre	Carbon Fibre	Carbon Fibre
WARRANTY	FRAME: 5-YEAR ELECTRONIC COMPONENTS: 1-YEAR BATTERY: 6-MONTH	FRAME: 5-YEAR ELECTRONIC COMPONENTS: 1-YEAR BATTERY: 1-YEAR	FRAME: 5-YEAR ELECTRONIC COMPONENTS: 1-YEAR BATTERY: 1-YEAR
COLOUR/FINISH	Gloss Carbon	Gloss Carbon & Graphite Gloss Carbon & Lilac Gloss Carbon & Teal	Gloss Carbon