

PORTABLE HYGIENE AID

PLEASE READ BEFORE USING YOUR PRODUCT

PRODUCT CODE: **BA02100006**

Single-use, portable aid to support independence and provide convenience in daily living.

DISPOSAL

- Ensure the aid is emptied fully after use.
- Maintain proper hygiene when handling the used aid.
- Dispose of the aid in the appropriate waste-disposal bin.

MEASUREMENTS (MM):	200 x 113 x 115
MAXIMUM CAPACITY (ML):	900
MATERIAL:	Paper Pulp

Due to manufacturing tolerances and continued product improvement, this specification can be subject to a variance of (+ or -) 3%.



WARNING

- Single use only.
- Do not exceed the maximum capacity.
- Take care when transporting the aid to avoid spills and structural damage.
- Intended for urine only – not suitable for other waste.